

# We protect you, so you can protect us.



***Introducing a revolution in respiratory protection—the most breathable, long-lasting N95 respirator ever designed.***

ElastoMaskPro is NIOSH-approved. It combines the durability of an elastomeric respirator with the simplicity of a disposable N95 mask. Soft, lightweight silicone is integrated with Forever Fit technology to provide all-day comfort, a secure seal, and proper fit even after 100,000 simulated uses.

A major technological ***and cost-saving advantage*** of the ElastoMaskPro is Filter Puck technology. Filter Pucks are easier to breathe through than what you're wearing now. One ElastoMaskPro is equal to thousands of disposable N95 filtering facepiece respirators.

ElastoMaskPro protects you, so you can protect us.

- ***Reduces waste—equivalent to hundreds of disposables for max savings.***
- ***Maintains fit even after 100,000 simulated uses.***
- ***Cleans easily with common disinfectants.***
- ***Provides all-day comfort and a secured seal (eyeglasses don't fog).***
- ***3 to 5 times easier to breathe through than many elastomeric respirators and 2 to 3 times more than many disposable N95 respirators.***
- ***Unique serial number and bar code ensures each user can maintain control of their respirator.***
- ***Low-profile design for improved field of vision.***




## ElastoMaskPro

NIOSH-approved

[Reusable-Respirators.com](https://www.Reusable-Respirators.com)

ElastoMaskPro™, Filter Puck™, Forever Fit™ are trademarks of Reusable Respirators, LLC.  
©2022 Reusable Respirators, LLC

# Advanced simplicity.



1. Shields protect filters and can be easily wiped down.
2. Extended vents allow for air flow and provide an easy user seal check.
3. Each respirator is printed with a unique serial number and barcode to allow for identification and disinfection tracking.
4. Each respirator can be cleaned with a simple wipe.
5. Forever Fit technology features integrated straps, allowing for thousands of uses with no degradation in fit performance.
6. Facepiece is available in two sizes and compatible with most users.
7. Replaceable Filter Puck technology provides the most breathable respirator on the market.

## Compare ElastoMaskPro...

Feature / Benefit	ElastoMaskPro	Traditional N95 FFR	Traditional Elastomeric Half-Mask Respirator
Reusable	✓	✗	✓
Minimum N95 protection	✓	✓	✓
Filters user's breath - no exhalation valve	✓	✓	✗
Simple, intuitive design	✓	✓	✗
Smooth surfaces to easily clean and disinfect (including straps)	✓	✗	✗
Easy user seal check without cross-contamination	✓	✗	✓
Lightweight	✓	✓	✗
Easy to breathe through	✓	✗	✗
Compact design to minimize visual impairment	✓	✓	✗
Designed for frontline workers	✓	✗	✗
Designed to meet <a href="https://www.cdc.gov/niosh/nppt/hospresptoolkit/pdfs/ProjectBREATHE-final-report-508.pdf">Project B.R.E.A.T.H.E. recommendations</a> <sup>1</sup>	✓	✗	✗

<sup>1</sup> <https://www.cdc.gov/niosh/nppt/hospresptoolkit/pdfs/ProjectBREATHE-final-report-508.pdf>

Scan to see more about ElastoMaskPro



[Reusable-Respirators.com](https://Reusable-Respirators.com)

ElastoMaskPro™, Filter Puck™, Forever Fit™ are trademarks of Reusable Respirators, LLC.  
©2022 Reusable Respirators, LLC